



High Performance Laminates for Reusable Gowns and Drape Systems

Stedair® TX L3

Stedair® isolation barriers specifically designed for protection against liquid penetration while remaining breathable to provide comfort for First Responders. TYPE 1 - Polyester with ePTFE film.

Product Benefits

Resistant to penetration of liquids


- *Wash durable - 75 wash sterilization cycles

- High breathability

SPECIFICATION	ANSI PB70:2012	RESULTS
Specific Weight		3.8 oz/yd²
CAN/CGSB-4.2 Nº 5.1		
Tearing Strength		W: 30 lbf F: 20 lbf
ASTM D1424		
Tensile Strength		W: 70 lbf F: 70 lbf
ASTM D5034		
Water Absorption	≤1.0 g	PASS
AATCC 42		
Water Resistance	10 min 10 psi	PASS
CAN/CGSB-4.2 Nº 26.5		
Water Vapor Resistance		5 000 g/m²/24 hrs
ASTM E-96 method BW		
Hydrostatic Pressure	≥50cm	PASS
AATCC 127		

WASHING RECOMMENDATIONS

- Verify all the equipment to be sure that there is no sharp object which could damage the fabric.
- Do not overload the equipment; 75-80% of the capacity is optimal.
- Temperature should not exceed 140°F for washing operations.
- Do not use Sodium Hypochlorite (beaching agent) as a disinfecting agent, use Hydrogen Peroxide (35%) in a concentration of 150 ppm, PH 10.5 maximum.
- High pressure extractor unit of type single stage plunger can damage fabrics made with Stedair® TX L3.
 Avoid excessive high temperature, or too much drying during the process. Make sure that you have an adequate cooling period and fold goods immediately after finishing drying to avoid wrinkling.
 Sterilization temperature should not exceed 270°F. It is highly recommended to implement a Quality
- Insurance program to ensure ongoing reliability in the performance of product containing Stedair® TX L3.



The information set forth herein reflects the performance of Stedair® TX L3 under controlled laboratory conditions.